



TYPE III
AMBULANCE

DEMERS TYPE III AMBULANCE

In saving lives, there is no such thing as compromise. There is no "good enough." You go all out. And so do we. Our ambulances meet or exceed North American standards for quality and safety, and come loaded with unique and innovative features. Demers is your trusted partner on the road – tested, hardworking and ready to roll.

Our Type III Ambulances are built for the long run. More than just an ambulance, our Type III vehicles are a top value for your investment, leading the industry for total cost of ownership and best resale value.

Built for both crew comfort and occupant safety, Demers Ambulances provide innovative standard features to make your job easier. With superior overall design and craftsmanship, Demers Type III vehicles combine safety, innovation and quality – all with the goal of getting you home safely.



"Demers sticks out as being the top of the line for safety, construction, and design."

- Brandon Mitchell, Administrative and financial chief for Durham County Paramedics (NC)

TYPE III ADVANTAGES



Affordably priced & loaded with features –
good for your bottom line



Top of the category in total cost of ownership



Designed for efficient patient care



Up to 8% fuel cost savings



Significant reduction in maintenance costs



Optional Mobility Track Seating

for better access to critical equipment



High-visibility flush mount windows




Standard multiplex electronic management system

for remote diagnostics, operational efficiency and reduction of down time



REQUEST A DEMO TODAY AT DEMERS-AMBULANCES.COM



For complete vehicle information including specs, drawings and standard options, visit demers-ambulances.com

TYPE III **SPECS**

	MX 151	MX 152	MX 164	MX 170
Overall Length	259" (657 cm)	264" (670 cm)	271" - 278" (688 cm - 706 cm)	278" -284" (706 cm - 721 cm)
Overall Width*	88" (223 cm)	86" (218 cm)	95" (241 cm)	95" (241 cm)
Overall Height**	103" (262 cm)	108" (274 cm)	106" (269 cm)	106" (269 cm)
Headroom	68" (173 cm)	72" (183 cm)	72" (183 cm)	72" (183 cm)
Wheelbase	138" (351 cm)	144" (366 cm)	158" -159" (401 cm - 404 cm)	158"-159" (401 cm - 404 cm)
Chassis	Ford E350	Mercedes Benz 3500	Chevrolet Express 3500/4500 Ford E350/450	Chevrolet Express 4500 Ford E450
Payload Capacity	Up to 2500 lbs. (1132 kg)	Up to 2300 lbs. (1043 kg)	Up to 4500 lbs (2041 kg)	Up to 4000 lbs. (1814 kg)

*Total width does not include the fenders.

**Total height of vehicle may differ according to chassis and suspension.



DEMERS MX 151 AMBULANCE

Also referred to as our “door forward” design, Demers MX 151 Ambulance features a curbside forward entry door design with the ALS equipment moved closer to the paramedics, keeping everything within arm's reach and enhancing patient care. Offering more standard features than other vehicles in its class, Demers MX 151 Ambulance is an affordably priced option.

- Door forward design
- Affordably priced and loaded with features
- Shorter wheelbase for urban response and maneuverability



“Demers ability to combine safety with customer needs is what puts them at the top. We have a partnership in the Build. It's the new Standard!”

- Fred Savage, Manager at Cox Health Systems

DEMERS MX 152 AMBULANCE

With distinctive details including a unique aerodynamic front cab riser and rear spoiler perfectly engineered to match the Mercedes Sprinter design, Demers MX 152 Ambulance offers eye-catching style. Offering our exclusive fold out entry step and a roomy interior with expanded storage compartments, the MX 152 combines dynamic looks with exemplary performance.

- Unmatched fuel economy
- Exclusive fold out entry step
- Expanded storage compartments





DEMERS MX 170 AMBULANCE

Demers MX 170 is a full-size Type III ambulance, offering a roomy interior for extensive medical equipment and an optimal ergonomic work environment for EMS specialists. Tested to the harshest global standards, this ambulance is built with unique features including full length, all-aluminum extruded door construction for maximum storage access without compromising structural integrity. Equipped with Demers mobility track seat, the MX170 is truly a category leader.

- > Loaded with standard features
- > Expanded ALS cabinet with ample storage
- > Cab extensions for crew comfort

DEMERS MX 164 AMBULANCE

The MX 164 is Demers workhorse and overall best value Type III ambulance. With thousands of units in service around the world, this design offers best in class safety and ergonomics, reduced maintenance and multiple choices of optional layout features. Offering our industry-leading mobility seat as an option, Demers MX 164 ambulance is one-of-a-kind.

- > Rugged
- > Reliable
- > Durable
- > Cab extensions for crew comfort layouts



“How Demers builds ambulances is unique to this industry. We found that the repeatability and the quality was outstanding.”

- Eric Scheumann, Las Vegas Fire & Rescue

THE DEMERS DIFFERENCE

Never compromise. It's our commitment to uphold unmatched industry standards for quality and safety. Our designs are more innovative. Our module is stronger. Our testing is more rigorous. Our ambulances are expertly crafted and designed for efficiency, providing the lowest total cost of ownership in the industry. That's a big deal, especially when you consider all that our vehicles have to offer. So much rides on our reputation. But what matters most and always has since 1960 is the safety of those who ride within.



55 YEARS
WORLD CLASS
AMBULANCES

**REQUEST A DEMO TODAY AT
DEMERS-AMBULANCES.COM**

28 Richelieu,
Beloeil, QC, J3G 4N5, Canada
info@demers-ambulances.com

demers-ambulances.com
(800) 363-7591

