



DEMERS MX 170

Type III Ambulance

In saving lives, there is no such thing as compromise. There is no “good enough.” You go all out. And so do we. Our ambulances meet or exceed North American standards for quality and safety, and come loaded with unique and innovative features. Demers is your trusted partner on the road – tested, hardworking, and ready to respond.

Our Type III ambulances are built for both crew comfort and occupant safety, providing innovative standard features to make your job easier. With superior overall design, Demers Type III vehicles combine safety, innovation, and quality – all with the goal of getting you home safely.

Demers MX 170 is a full-size Type III ambulance, offering a roomy interior for extensive medical equipment and an optimal ergonomic work environment. Tested to the harshest global standards, this ambulance is built with unique features including all-aluminum extruded door construction creating maximum storage access, and can be equipped with Demers mobility track seat amplifying crew/patient safety.

- Designed for efficient patient care
- Cab extended into module for more cab space & improved ergonomics
- Expanded ALS cabinet with ample storage
- Dual HVAC ducted system
- Aerodynamic roof design with incorporated full LED warning lights for increased intersection visibility (optional)



Demers-Ambulances.com

info@demers-ambulances.com

(800) 363-7591





DEMERS MX 170 FEATURES

-  **Designed for efficient patient care** with customizable innovative options to suit your department's needs
-  **Optional Mobility Track Seating** for better access to critical equipment while seated and belted
-  **Exclusive HVAC airflow system** with separate hot and cold air ducting for improved patient comfort
-  **Aluminum clever cabinetry** is pull tested for occupant safety and interior durability
-  **Top of the category in total cost of ownership** built with lightweight materials and ECOSmart features
-  **High-visibility flush mount windows** prevents leaks and lowers repair costs
-  **Standard multiplex electronic system** for remote diagnostics, operational efficiency and reduction of down time

DEMERS MX 170 SPECS

Chassis	Chevrolet Express 4500 Ford E450
Drive Type	2x4
Overall Length	278" to 284"
Overall Width*	95"
Overall Height**	106"
Interior Headroom	72"
Wheelbase	158" to 159"
Avg. Payload Capacity	Up to 4,000 lbs.

*Total width does not include the fenders.

**Total height of vehicle may differ according to chassis and suspension.



ECOSMART FEATURES

Economically and Ecologically Smart

Demers ECOSmart Technologies reduce operating costs, fuel consumption and impact on the environment.

Aerodynamic roof design, up to 8% savings in fuel costs

Demers Ambulances can come equipped with the signature aerodynamic roof design with incorporated LED warning lights for enhanced safety and durability, saving fuel costs.

Lightweight materials, for better fuel economy & higher resale value

Every detail counts when it comes to improving your performance and driving more money to your bottom line. Demers Ambulances are engineered to be as light and durable as possible so that you can carry the equipment you need while lowering fuel costs.

Anti-idling system: reduces your operational costs, improves environment imprint

Powered by the Demers Electronic Management System, the Demers ECOSmart system creates an intelligent, anti-idling system for the vehicle's engine—an innovation that can deliver fuel savings. This also translates into a reduction of CO₂ emissions.

THE DEMERS DIFFERENCE

Never compromise. It's our commitment to uphold unmatched industry standards for quality and safety. Innovative designs, strong modules, and rigorous testing. Our ambulances are expertly crafted and designed for efficiency, aiming to provide the lowest total cost of ownership in the industry. So much rides on our reputation. But what matters most, and always has since 1960, is the safety of those who ride within.



Brought to you exclusively by
Demers Ambulances USA Inc. and
their authorized Dealer partners.

Demers-Ambulances.com

info@demers-ambulances.com

(800) 363-7591

