



## DEMERS MXP 150

### Type I Ambulance

When saving lives, there is no such thing as compromise. There is no 'good enough'. You go all out. And so do we. Demers Type I Ambulances are built for the long haul, handling extreme and demanding conditions – tested, hardworking, and ready to respond.

The Demers MXP 150 Ambulance is built to handle a wide range of conditions and environments with its rugged exterior and top-of-the-line design details, including an industry-leading insulation package. It blends ergonomics, innovative features & functionality to elevate the overall patient care experience with standard dual ducted HVAC, aluminum clever cabinetry, and optional mobility track seating.

More than just an ambulance, our vehicles are a top value for your investment, a leader in the industry for total cost of ownership and best resale value. Demers is your trusted partner on the road – safety, innovation, and quality – all with the goal of getting you home safely.

- Spacious front cabin with pass-through window
- Spacious compartment (ALS) on front wall with adjustable shelves
- Dual HVAC ducted system
- More standard features than comparable Type I Ambulances
- Aerodynamic roof design with incorporated full LED warning lights for increased intersection visibility (optional)







MXPI50-TI-04-2025

[Demers-Ambulances.com](http://Demers-Ambulances.com)  
(800) 363-7591

[info@demers-ambulances.com](mailto:info@demers-ambulances.com)  
f X in @ YouTube



## DEMERS MXP 150 FEATURES

-  **More standard features than comparable Type I ambulances**
-  **Engineered for extreme conditions**  
from desert heat to frigid environments
-  **Exclusive HVAC airflow system**  
with separate hot and cold air ducting for improved patient comfort
-  **Numerous customizable interior options**  
allowing you to create an ambulance that meets your service needs
-  **Best-in-class payload capacity**  
up to 4,500 lbs. giving you room for upgrades and storage weight
-  **High-visibility flush mount windows**  
prevents leaks and lowers repair costs
-  **Full LED warning light packages**  
to take on the darkest nights and brilliantly light the way

## DEMERS MXP 150 SPECS

Chassis	Dodge 4500   Ford F350/450/550
Drive Type	2x4 or 4x4
Overall Length*	285" to 291"
Overall Width	95"
Overall Height*	112"
Interior Headroom	72"
Wheelbase*	168.5" to 171"
Avg. Payload Capacity	Up to 4,500 lbs.

\*Depending on chassis selection, total height of vehicle may differ according to chassis and suspension.



## ARCTIC RANGE HVAC

**Arctic Range is an optional premium insulation package & climate control system**

- Industry's most efficient insulation package
- Extreme Climate Tested
- Reaching lowest interior temperatures at exceptional speed

## ECOSMART FEATURES

Economically and Ecologically Smart

**Demers ECOSmart Technologies reduce operating costs, fuel consumption and impact on the environment.**

**Aerodynamic roof design, up to 8% savings in fuel costs**

Demers Ambulances can come equipped with the signature aerodynamic roof design with incorporated LED warning lights for enhanced safety and durability, saving fuel costs.

**Lightweight materials, for better fuel economy & industry leading operating costs**

Every detail counts when it comes to improving your performance and driving more money to your bottom line. Demers Ambulances are engineered to be as light and durable as possible so that you can carry the equipment you need while lowering fuel and maintenance costs.

**Anti-idling system: reduces your operational costs, improves environment imprint**

Powered by the Demers Electronic Management System, the Demers ECOSmart system creates an intelligent, anti-idling system for the vehicle's engine—an innovation that can deliver fuel savings. This also translates into a reduction of CO<sub>2</sub> emissions.

## THE DEMERS DIFFERENCE

Never compromise. It's our commitment to uphold unmatched industry standards for quality and safety. Innovative designs, strong modules, and rigorous testing. Our ambulances are expertly crafted and designed for efficiency, aiming to provide the lowest total cost of ownership in the industry. So much rides on our reputation. But what matters most, and always has since 1960, is the safety of those who ride within.



Brought to you exclusively by  
Demers Ambulances USA Inc. and  
their authorized Dealer partners.

[Demers-Ambulances.com](http://Demers-Ambulances.com)

**(800) 363-7591**

[info@demers-ambulances.com](mailto:info@demers-ambulances.com)

