

FXP 174

Type I Ambulance



WHERE INNOVATION MEETS SUPERIOR PATIENT CARE

Experience the future of emergency response with a groundbreaking Type I ambulance that redefines crew safety, performance, and patient care.

Purpose-built for first responders, it features an ergonomic design informed by real EMS field experience and ambulance dealer insights.

The state-of-the-art clinical layout and exclusive FX Action Seats enable superior patient care, when it matters most.



PRODUCT OVERVIEW

- › The revolutionary FXP 174 ambulance transforms emergency response with significant advancements & improvements, maximizing crew safety and patient outcomes
- › The Type I model pushes the boundaries of conventional ambulance layouts while integrating innovative technologies
- › Purpose-built for first responders, this ambulance optimizes crew safety and comfort through ergonomic design shaped by extensive industry feedback
- › The patient compartment features state-of-the-art clinical design with ergonomics that optimize crew workflow and mobility while providing unmatched access to patients
- › Exclusive FX Action Seats ensure EMS safety while enabling optimal mobility and patient care

PRODUCT VALUES

- › **Provider Wellbeing:** Enhanced safety, improved comforts and care for EMS crew and patients
- › **Operational Excellence:** Optimized workflows enable faster, safer emergency response
- › **Safety Integration:** FX Action Seats promote crew protection without compromising vital patient care
- › **Modern Technology:** Delivers cutting-edge solutions that attract next-generation EMS professionals
- › **Innovative Configuration:** Best-in-class layout maximizes space & functionality

INTENTIONALLY DESIGNED AND THOUGHTFULLY ENGINEERED FOR YOU.



FXP174-FACT-09-2025

Demers-Ambulances.com
(800) 363-7591

info@demers-ambulances.com



A REIMAGINED AMBULANCE

FEATURING A NEVER-BEFORE-SEEN CLINICAL DESIGN

FX ACTION SEATS

Revolutionizing EMS safety & patient care with proprietary technology

- Purpose-built ergonomic seat design for First Responders
- Maximizes mobility while enabling optimal patient care with the push of a button
- More opportunities to remain safely buckled during patient care
- Stabilized seat to help perform critical tasks
- Intuitive design enables Medics to maneuver using legs, releasing hands to care for patients

- Installed near stretcher & equipment to prevent unbuckling
- Enhanced safety with emergency lock (26G acceleration)
- Belt & occupant sensor for driver awareness



Revolutionary new FX Action Seats offer medics freedom of movement and accessibility to patient & equipment while belted.

PRODUCT FEATURES

Patient Compartment Specifications

- 174" Length x 95" Width
- 72" Interior Headroom
- Door Forward Design

Chassis Configurations

- Ford F450 & F550, Gas & Diesel

Vehicle Exterior

- Door Forward Design enhances ergonomics, workflow and space efficiency
- Taller rear doors enables efficient entry & exit
- Advanced lighting system with breakthrough flash patterns provides industry-leading visibility
- Windows on both curbside & streetside enhance crew visibility
- Double-door compartment maximizes storage capacity
- Triple D-bottle storage saves space & accessible via C2 curbside door
- Exterior compartment with convenient access for electrical system maintenance

Patient Compartment (Module)

- Clinical Layout prioritizes ergonomics, safety and comfort with state-of-the-art design
- Functional Interior Storage with improved access to interior compartments and workstations, optimized to keep paramedics seated while accessing patient and equipment
- One Click 4-Point Seatbelt allows freedom to stand while remaining buckled
- Ultrawide Electric Passthru Window optimizes visibility & reduces noise with sealed barrier
- Linear lighting with adjustable intensity and colors enhancing patient care & comfort
- Dual HVAC Ducted Diffusion System provides separate hot and cold air ducting for improved patient & crew comfort

CONTACT INFORMATION

Kris Turner, Director Business Channel 419-203-8639
Vanessa Bomboir, Director Marketing 306-221-6584

